

INVENT – A DECISION SUPPORT SYSTEM FOR MANAGING INSPIRED OXYGEN – PROSPECTIVE EVALUATION IN AN INTENSIVE CARE UNIT

Karbing, DS<sup>1</sup>; Allerød, C<sup>1,2</sup>; Thorgaard, P<sup>2</sup>; Andreassen, S<sup>1</sup>; Kjærgaard, S<sup>2</sup> & Rees, SE<sup>1</sup>

<sup>1</sup>Center for Model-based Medical Decision Support, Department of Health Science and Technology, Aalborg University, Aalborg, Denmark; <sup>2</sup>Anaesthesia and Intensive Care, Region North Jutland, Aalborg Hospital, Aarhus University, Aalborg, Denmark.

Introduction: INVENT is a model-based decision support system for ventilator management providing suggestions of tidal volume, respiratory frequency and FiO<sub>2</sub> [1]. As a first step, this study prospectively evaluated INVENT for managing FiO<sub>2</sub> in an ICU.

Methods: 13 patients were available for analysis with 45 experiments where both attending clinicians and INVENT managed FiO<sub>2</sub>. Clinicians set FiO<sub>2</sub> as in normal practice. INVENT estimated model parameters and then suggested the FiO<sub>2</sub> representing the optimal compromise between oxygenation goals and minimizing FiO<sub>2</sub> to avoid adverse effects of hyperoxia [2, 3].

Results: Median (range) INVENT and clinician changes in FiO<sub>2</sub> from baseline were -0.03 (-0.16 to 0.12) and 0.00 (-0.10 to 0.05). Median SaO<sub>2</sub> changes from baseline for INVENT and clinicians were 0.00 (-0.04 to 0.05) and 0.00 (-0.05 to 0.03). INVENT FiO<sub>2</sub> and SaO<sub>2</sub> ranges were 0.26–0.54 and 0.94–0.99. Clinician FiO<sub>2</sub> and SaO<sub>2</sub> ranges were 0.25–0.70 and 0.91–0.99. Intra-patient SaO<sub>2</sub> variation across experiments was at most 0.02 in 12 patients for INVENT compared to 9 for clinicians.

Discussion: INVENT SaO<sub>2</sub> was within previously defined safe limits for ALI/ARDS [4]. INVENT accomplished narrower FiO<sub>2</sub> and SaO<sub>2</sub> distributions than clinicians, i.e. standardizing FiO<sub>2</sub> management, and better maintained adequate SaO<sub>2</sub> when patients changed status. INVENT may enable automatic regular reevaluations of patients improving care and freeing clinicians' focus.

REFERENCES

1. Rees SE et al. Using physiological models and decision theory for selecting appropriate ventilator settings. *J Clin Monit Comput* 2006; 20:421–9.
2. Altmeier WA, Sinclair SE. Hyperoxia in the intensive care unit: why more is not always better. *Curr Opin Crit Care* 2007; 13:73–8.
3. Aboab J et al. Effect of inspired oxygen fraction on alveolar derecruitment in acute respiratory distress syndrome. *Intensive Care Med* 2006; 32:1979–86.
4. The ARDS Network: Ventilation with lower tidal volumes as compared with traditional tidal volumes for acute lung injury and the acute respiratory distress syndrome. *N Engl J Med* 2000; 342:1301–8.